


VAROS 307 Firmware-Update

Before starting the update be sure that the instrument is connected to the power supply. If you complete the update with your laptop then this should also be connected to the power supply.


1. Open the workstation




2. Double-click for e.g. CD drive C




3. Using the right-click make a new folder – new – folder




4. Enter a name (e.g. Update firmware 307)




5. Save the ZIP data in the new folder




6. Double-click the ZIP file




7. Mark all the data and then choose copy (edit – copy)




8. Click the return button and then choose insert (edit – insert)



9. Double-click the update software (exe.data)




10. Select the serial connection of your PC which is connected with the VAROS 307.




Attention:

For this connection you need a serial extension cord (9pol SUB-D) with a null-modem allocation or a USB serial adaptor with a fit on null-modem adaptor.


11. Click the button „ Update Firmware ”




12. In the dialog for Firmware-data choose the BIN-data (e.g. Vxx_04d.bin) by double-clicking.



13. The update will start automatically after choosing the Firmware-data.



14. Please wait until the update has been successfully installed.



IF THE CONNECTION IS INTERRUPTED BEFORE THEN IT MIGHT DAMAGE THE INSTRUMENT: IN THIS CASE IT WOULD HAVE TO BE SENT TO US AT THE OWNERS EXPENSE: